

Julabo Case Study

JULABO PRESTO® W92tt

Cooling and heating a 100 liters reactor
between -50 °C and +100 °C



Objective

This case study tests the reproducibility of heat-up and cool-down processes of a JULABO PRESTO® W92tt with a 100 liters glass reactor. The W92tt is connected to the reactor via two 2.0 m metal tubings. The W92tt is programmed to cycle between -50 °C and +100 °C.

Test Conditions

JULABO unit	JULABO PRESTO® W92tt
Cooling power	+20 °C 19.0 kW
	0 °C 15.5 kW
	-20 °C 9.5 kW
Heating capacity	36 kW
Band limit	70 K
Flow pressure	0.33 bar
Bath fluid	JULABO Thermal HL80
Reactor	100 liters glass reactor (Büchiglas) filled with 100 liters Thermal HL80
Control	External (ICC)

Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	3 x 400 V / 50 Hz



Test Results

See chart on back page: The reactor was repeatedly heated up and cooled down between -50 °C and +100 °C with a PRESTO® W92tt. Heat-up and cool-down processes were repeated exactly several times.

Tip


You can also use the robust Pt100 with PTFE coating.

More tips on
back page >>



JULABO GmbH
Eisenbahnstraße 45
77960 Seelbach / Germany
Tel. +49 (0) 7823 51-0

Measured with EasyTEMP Professional



Tip

Make use of the option to regulate the pump pressure. You can define the desired pressure in the PRESTO® settings.



Tip

The Ethernet interface permits full access to all operational functions of the PRESTO®.



JULABO GmbH
 Eisenbahnstraße 45
 77960 Seelbach / Germany
 Tel. +49 (0) 7823 51-0